

DigiDop Dopplers

"Built to Last"

TOUGH and *Sensitive*

Minimize Service Returns

Flexible, magnetic probe holder won't crack or break. Tough stainless/Polycarbonate case.

Configure to your application -

Interchangeable vascular & obstetrical probes. Smooth. Comfortable. Easy clean.

Five Year Warranty- Includes Probes

Zero Cost Warranty includes parts and labor on Doppler and probes.

Sensitive - More listening, less searching

Unique digitally Optimized Processing (Dop™) provides up to 6x signal of other Dopplers.

Save \$100's in service costs

User replaceable cords, great warranty



MADE IN

THE USA



What makes the DigiDop difference!

① Best Sound – Digitally Optimized Processing

The DigiDop provides a unique powerful combination of tough durability along with proven, measured sensitivity. We know that the real key to a great Doppler is performance – sensitivity and sound quality. While others may say this; we measure it – and have found that our signal is more than six times greater than some other Dopplers. Our Digitally Optimized Processing - DOP™ - eliminates most background noise while letting the Doppler sounds through bright and clear. With DOP each interchangeable probe can be optimized with different filtering in the main unit.

② Best Warranty -5 year \$0 cost warranty -ALL parts. ALL labor.

True Lower Costs of Ownership – Five Year \$0 Warranty – Including Probes Our 5 year warranty covers parts and labor, including our probes and cords. This can easily be a savings of hundreds of dollars compared to other Dopplers on a single service issue. We are able to offer this warranty because our design and tough materials significantly reduce the number and cost of warranty issues.

③ Best Support & Service – Knowledgeable Staff & Guaranteed Service

Description	Model	5 Year \$0 Warranty	Unbreakable Probe Holder	Interchangeable Probes?	Display?	Rechargeable?
DigiDop 300	DD-300-**	✓	✓	✓		
DigiDop 301	DD-301-**	✓	✓	✓		✓
DigiDop 700	DD-700-**	✓	✓	✓	✓	
DigiDop 701	DD-701-**	✓	✓	✓	✓	✓

** When ordering, place the desired probe number from the probe section below in this location.

Probe Options

D2	2MHz obstetrical probe - for deeper penetration
D2W	2MHz waterproof obstetrical probe for deeper penetration/water applications
D3	3MHz obstetrical probe - early heartbeat detection & general use
D3W	3MHz waterproof obstetrical probe for general use water applications
D5	5MHz vascular probe - pen style vascular for deeper vessels
D8	8MHz vascular probe - pen style for superficial vessels & general use

Technical Specifications

Weight	340 grams (12 ounces)
Dimensions	150 x 65 x 35 mm (6.0 x 2.5 x 1.25 in)
Battery Life	900 minutes Alkaline/720 minutes NiMH
Warranty	5 years ALL parts ALL labor

Obstetrical Probes – Sensitive, Durable, Versatile

Our interchangeable obstetrical probes are easy to hold and clean. Each Doppler probe is individually tuned for maximum sensitivity. The probe cords are user replaceable.



D2

2MHz – Deeper penetration



D2W

Waterproof 2MHz for water bath/birth



D3

3MHz – Very sensitive to early fetal heart, best all around OB probe



D3W

Waterproof 3MHz probe for water bath/birth

Vascular probes are also available:

D5 – 5MHz, D8 – 8MHz, PPG – audio ppg

DigiDop Dock Options

Secure your DigiDop by mounting on the wall, an IV pole, a roll-stand, or a table dock. Theft resistant fastener keeps it there.



Roll Stand Option
includes basket



Table Dock Option



IV Pole Option



Wall Option

How to order!

Simply add probe choice to end of main unit choice
For example, DD-700-D3 is a Display Doppler with 3MHz probe